

## **Supplementary information**

# **Restraint Stress-Induced Immunosuppression Is Associated with Concurrent Macrophage Pyroptosis Cell Death in Mice**

**Chi-Cheng Li<sup>1,2</sup>, Rina Munalisa<sup>3</sup>, Hsuan-Yun Lee<sup>3</sup>, Te-Sheng Lien<sup>3</sup>, Hao Chan<sup>3</sup>, Shih-Che Hung<sup>3</sup>, Der-Shan Sun<sup>3</sup>, Ching-Feng Cheng<sup>4,5</sup> and Hsin-Hou Chang<sup>3,\*</sup>**

<sup>1</sup> Department of Hematology and Oncology, Buddhist Tzu Chi General Hospital, Hualien 970, Taiwan; kevinlcc1234@gmail.com

<sup>2</sup> Center of Stem Cell & Precision Medicine, Hualien Tzu Chi Hospital, Hualien 970, Taiwan

<sup>3</sup> Department of Molecular Biology and Human Genetics, Tzu-Chi University, Hualien 970, Taiwan; 108727110@gms.tcu.edu.tw (R.M.); 109727102@gms.tcu.edu.tw (H.-Y.L.); alan211@mail.tcu.edu.tw (T.-S.L.); chanhao1011@gmail.com (H.C.); 102353113@gms.tcu.edu.tw (S.-C.H.); dssun@mail.tcu.edu.tw (D.-S.S.)

<sup>4</sup> Department of Pediatrics, Taipei Tzu Chi Hospital, Buddhist Tzu Chi Medical Foundation, Taipei 231, Taiwan; chengcf@mail.tcu.edu.tw


<sup>5</sup> Institute of Biomedical Sciences, Academia Sinica, Taipei 115, Taiwan

\* Correspondence: hhchang@mail.tcu.edu.tw; Tel.: +886-3-8565301 (ext. 2667)

**Key words:** Restraint stress, immunosuppression, macrophage, cell death, pyroptosis, psychological stress, ambient cold exposure, intravenous immunoglobulin

**Figures S1:** pages 2

Supplemental Figure S1



**Figure S1.** Restraint stress suppressed the live phagocyte function as indicated by the suppression on the engagement between live splenic macrophage and fluorescent nanoparticles (NPs) in mice. (A) Experiment outline, and (B) percentages of NP-engaged live cells are indicated. Live splenocyte populations of respective groups were normalized to 100 %. \*  $P < 0.05$ , vs. respective 0 h groups (untreated); #  $P < 0.05$ , vs. respective no stress (−restraint stress) groups. n = 6 (three experiments with total 6 mice per group). Live cell populations were determining using a live/dead cell labeling kit (Zombie NIR™ Fixable Viability Kit, Biolegend).