

Supplemental material to:

Treatment of early allergic and late inflammatory symptoms of allergic rhinitis with Petasites hybridus leaf extract (Ze 339): Results of a non-interventional observational study in Switzerland




Figure S1: Gas chromatogram of *Petasites hybridus* leaf extract Ze 339 (batch 091376). Quantitative determination of total petasins (petasin, isopetasin, neopetasin), the active compounds of Ze 339 using gas chromatography and a flame ionization detector (FID). GC-column 100% polydimethylsiloxane (e.g. DB-1, length: 25 m, ID: 0.32 mm, dF: 0.52 μ m); Injector temperature: 270 °C; Injection volume 1 μ L.