



## Supplementary materials

**Table S1.** The sequences and annotations of promoters used in this study.

Promoter	Sequence	Source
cbp_2	GAGTCGTGACTAAGAACGTCAAAGTAATTAAACAATACAGCTATTTCTCATGCTTT ACCCCTTCATAAAATTAAATTATCGTTATCATAAAAAATTATAAGCGTTATATTG CTTGCCTGGATATAGTGCTGGCATTGTTGGCAGGAAATGTCGGAGTAAGGTGGA TATTGATTTGCATGTTGATCTATTGCAATTGAAATGATTAGTTATCCGTAATATTAAATT AATCATATCATAAAATTATATCATAATTGTTGACGAATGAAGGTTTGGATA AATTATCAAGTAAAGAACGCTAAAATTITGGCTAAAATATCAAATGACCACT <b>TGAATTAAATATGGTAAAGTAGATATAATATTGGTAAACATGCCCTCAGCAAGGTT</b> AGATTAGCTGTTCCGTAAATTACCGTATGGTAAACGGCAGTCAGAAAAATA AGTCATAAGATTCCGTATGAAATATACTTCGTAGTTAATAAGAGATATGAG GTAAGAGATACAAGATAAGAGATATAAGGTACGAATGTATAAGATGGTCTTTAG GCACACTAAATAAAAAACAAATAAACGAAATTAA <u>AGGAGGACGAAAG</u>	[1]
gapDH	TAATTACTGTATCTCTGGCATTGCCAGGTTAATAAAGATTAAATTATTGACTA GAAATAAAAAATTGCCATAATATTAAATGGACAAAAAAACAAAGAATTACATCA AAGGAAGATAAAAATACTTGTAAAAATTAAATTATTTTATCTAAACTATTGAA AATGAAAATAAAATAATAAAATGAATCATAGTCAAGAGATACTGCCAGAGGA TGAATATTACTGCATTATGGCAGCTAATAGAGGCATTAAATTAAATTAAATT TAATTACAAT <u>AGGAGGCGATATTAAATG</u>	[1]
gapDH/tetO	TAATTACTGTATCTCTGGCATTGCCActctatcattgatagag <u>TAAAATTATTGACTAGAA</u> ATAAAAAATTGCCATAATATTAAATGGACAAAAAAACAAAGAATTACATCAAAG GAAGATAAAAATACTTGTAAAAATTAAATTATTTTATCTAAACTATTGAA <u>actc</u> tatcattgatagag <u>TAAAATGAATCATAGTCAAGAGATACTGCCAGAGGATGAATATT</u> ACTGCATTATGGCAGCTAATAGAGGCATTAAATTAAATTAAATTACAA <u>TAGGAGGCGATATTAAATG</u>	[1], This paper
eno	GGAAATATTAAAATGAAATGTTGAAAAATGTTTAAGATGGCATTATGGATA AAATATACTATGGTTGCAATAATGCTTCTATTAAATTGGACTTGTGTAATATG GTAGAAGGATGCAGTGTAAATTAAACATATAAAATAAGCTATATGAAG <u>GGAG</u> AATGGAGA	[1]
0815	ATAATATATATCCGTATTTACAATATCCATACGAATATGGATTATTTATTTGTTAT ATTTTAACAAAAAAATTATTAAACCTTTCAATGGTGGTATTATATATATTGGTCA CAGTTCAATGAATC AAAAATAA <u>AGAGGTGTTAAT</u>	[1]
0966	TTCTTA <u>CTACTCCCTTGATCTACAGACAAATTTCCTCTTAACTTTCTGTAAACC</u> GTATCCGCATACTTTGTCAAAAGTAATTATAGATTATTACCGTCAAATCAAAT AATAAATTATATGCCAATTCTTGCATCCGTTAAAGTCTGTGCTATAATTATA TGGTAAAATATCATAAAAAAAT AAAATAAAAGTTAA <u>AGGA</u> ACTGATT	[1]
2638	GATAAACAAAGGACGGTCAGGGCTCTGCTCATCCTACTCTGCATTGAAAAAGGT AGGATGAATTTTATTAAATCTTATTGAAAAAAATTGAAACCGTTTATT <u>A</u> <u>AAAAAAAGTGGTATATTATAATAGCAATTGATTGGTAAAAAAATTAAAG</u> CAAACAGAATAATAACAAA <u>GTAAGGAGGA</u> ATTGTT	[1]
2926	AAAATATACAAAGGTTCTGTGTTTAATACCGTTATGTTAATATAATGTAATATAT ATTTTATAATAATATGTATGAGAGATAGTGTGCTATTGCTATAAGAAT <u>GAGG</u> <u>AGGAACTAG</u>	[1]
P <sub>xyl</sub> *	TTTGATATT <u>CCTCCTTAA</u> ATAATTGTAATACCTTACACAAAAATAAAGGTT ATTTCGATT <u>GACAAAGATAATTAAATTTATTAGTC</u> Ataagttagttaatataactaa	[2]

caaaaATAAAGCAAGTAAAATACCTAAATATAAAAAATTAGGATAGGAAAA  
CGATAGTTATGAAGTGGCATTCAAGGAGGGAT GCAT

Predicted features are marked as follows: ribosome binding sites are underlined, -35 boxes are bolded, -10 boxes are italicized and bolded, and operators tetO1 (in gapDH 2tetO1) and xylO (in P<sub>xyl</sub>) sequences are in lowercase.

\*P<sub>xyl</sub> contains two sets of -35/-10/RBS. The first RBS, -10, and -35 sequences are for the transcription of the *xylR* gene downstream of this set. The second RBS, -10, and -35 sequences are for the gene of interest, i.e., the gene whose transcription will be controlled by the promoter.

**Table S2.** List of primers used in this study. All primers were ordered from IDT.

Oligonucleotide	Sequence (5' - 3')	Purpose
adhB BglIII for	TAATCGAGATCTATGAAAGGTTTGAATGC	Amplification of <i>adhB</i> *
adhB Sall rev	TAATCGGTGACTTATGCTAATATTACAACAGGTTG	
Bgal BglIII for	TAATCGAGATCTATGAATGTGTTATCCTCAATTG	Amplification of <i>lacZ</i> *
Bgal Sall rev	TAATCGGTGACCTAACCTTCCGGCTTC	
PgapDH PspOMI for	TAATACGGGCCCTAATTACTGTATCTCTGGC	Amplification of P <sub>gapDH</sub> *
PgapDH BamHI rev	TAATACGGATCCTAATATCCCTCCTATTGTAAATTAAA	
P2638 PspOMI for	TAATCGGGGCCGATAAACAAAGGACGGTTC	Amplification of P <sub>2638</sub> *
P2638 BamHI rev	TAATCGGGATCCAACAAATTCCCTCCTACTTTG	
P2926 PspOMI for	TAATCGGGGCCAAATATACAAAGGTTCTTGTG	Amplification of P <sub>2926</sub> *
P2926 BamHI rev	ATAATCGGATCCAAATATACAAAGGTTCTTGTG	
Pcbp_2 PspOMI for	TAATCGGGGCCGAGTCGTGACTAAGAACG	Amplification of P <sub>cbp_2</sub> *
Pcbp_2 BamHI rev	TAATCGGGATCCCTTCGTCCCTTAAATTITTC	
Peno PspOMI for	TAATCGGGGCCCGAAATATTAAAATGGAAATGTTG	Amplification of P <sub>eno</sub> *
Peno BamHI rev	TAATCGGGATCCTCTCCATTCTCCCTCATATAG	
P0966 PspOMI for	TAATCGGGGCCCTTACTACTCCCTTGC	Amplification of P <sub>0966</sub> *
P0966 BamHI rev	TAATCGGGATCCAAATCAGTCCTTTAACTTTATT	
P0815 PspOMI for	TAATCGGGGCCATAATATATCCGTATTTACAATATCC	Amplification of P <sub>0815</sub> *
P0815 BamHI rev	TAATCGGGATCCATTAAACACCTCTTATTTTGATT	
PgapDH BamHI for	AATCAGGATCCGGCCCTAATTACTGTATCTCTGGC	Amplification of gapDH-lacZ module to add to pUC19
PgapDH BglIII for	AATCAAGATCTGGCCCTAATTACTGTATCTCTGGC	Amplification of P <sub>gapDH/tetO-lacZ</sub> module to add to pMTL86251 downstream of P <sub>eno-tetR</sub>
TetO1_1 for	CTATTGAAA <u>ACTCTATCATTGATAGAGTAAAATGAATC-ATAGTGC</u>	Add tetO1 to region between 2nd -35 and -10 in P <sub>gapDH</sub>
TetO1_1 rev	ACTCTATCAATGATAGAGTTCAATAGTTAGATAAAAAAT	
TetO1_2 for	TCCCTATCAGTGATAGATAAAATTATTGACTAGAAATAAAAAAATT	Add tetO1 to region between 1st -35 and -10 in P <sub>gapDH</sub>
TetO1_2 rev	GTCC TCAATAATTACTCTATCAATGATAGAGTTGGCAATGCCAGAGAG	
Eno rev EcoRI	ATAC	Amplification of P <sub>eno</sub> to clone in front of tetR
Eno for BamHI	ATTCTGAATTCTCTCCATTCTCCCTCATATAG	
tetR for EcoRI	ATAATCGGATCCGGAAATATTAAAATGGAAATGTTGAA	Amplification of tetR from
tetR rev PspOMI	ATCTGAATTCTAGATTAGATAAAAGTAAAGTG	E. coli to clone behind P <sub>eno</sub>
	ATTGGGCCCAACTCGACATCTGGTTAC	


\*Products were amplified from the appropriate pDGO plasmids [1] (Table S3).

**Table S3.** List of plasmids. All plasmids (with the exception of pDGO based plasmids) use pMTL86251 as the backbone.


Plasmid	Description	Source
pMTL86251	Modular shuttle vector for use in clostridial bacteria	[3]
pAL58	Plasmid containing adhB reporter gene driven by gapDH promoter	This paper
pAL61	Plasmid containing adhB reporter gene driven by 2926 promoter	This paper
pAL2B	Plasmid containing lacZ reporter gene driven by eno promoter	This paper
pAL60	Plasmid containing lacZ reporter gene driven by 2638 promoter	This paper
pAL62	Plasmid containing lacZ reporter gene driven by 2926 promoter	This paper
pAL64	Plasmid containing lacZ reporter gene driven by cbp_2 promoter	This paper
pAL66	Plasmid containing lacZ reporter gene driven by gapDH promoter	This paper
pAL731	Plasmid containing lacZ reporter gene driven by 0966 promoter	This paper
pAL751	Plasmid containing lacZ reporter gene driven by 0815 promoter	This paper
pAL88	Plasmid containing lacZ reporter gene driven by $P_{gapDH/tetO}$	This paper
pAL111	Plasmid containing lacZ driven by $P_{gapDH/tetO}$ and tetR driven by $P_{eno}$	This paper
pAL112	Plasmid containing lacZ driven by $P_{gapDH}$ and tetR driven by $P_{eno}$	This paper
pALX3	Plasmid containing lacZ driven by $P_{xyl}$ promoter and $xylR$ repressor gene	This paper
pDGO89	Plasmid containing adhB reporter gene driven by 0815 promoter	Gift from Dan Olson [1]
pDGO95	Plasmid containing lacZ reporter gene driven by gapDH promoter	Gift from Dan Olson [1]
pDGO102	Plasmid containing lacZ reporter gene driven by eno promoter	Gift from Dan Olson [1]
pDGO105	Plasmid containing lacZ reporter gene driven by 0966 promoter	Gift from Dan Olson [1]
pDGO106	Plasmid containing lacZ reporter gene driven by 2638 promoter	Gift from Dan Olson [1]
pDGO108	Plasmid containing lacZ reporter gene driven by 2926 promoter	Gift from Dan Olson [1]
pDGO117	Plasmid containing lacZ reporter gene driven by cbp_2 promoter	Gift from Dan Olson [1]

**Table S4.** Plasmids constructed for this study with promoter activity determined in *H. modesticaldum* and *E. coli*.

Plasmid	Promoter	Reporter gene	Detected activity in <i>H. modesticaldum</i> (U/mg protein)	Detected activity in <i>E. coli</i> (U/mg protein)
pAL2B	eno	lacZ	0.60 ± 0.019	0.17 ± 0.0050
pAL60	2638	lacZ	0.067 ± 0.0042	0.037 ± 0.0028
pAL62	2926	lacZ	0.39 ± 0.011	0.39 ± 0.011
pAL64	cbp_2	lacZ	0.64 ± 0.014	0.39 ± 0.011
pAL66	gapDH	lacZ	0.79 ± 0.027	0.25 ± 0.0069
pAL731	0966	lacZ	0.11 ± 0.0044	0.11 ± 0.0045
pAL751	0815	lacZ	0.10 ± 0.0044	0.013 ± 0.00077



**Figure S1.** The normalized β-gal activity in *H. modesticaldum*. The activity of β-gal in *H. modesticaldum* transformed with pALX3 and grown with xylose at the indicated concentrations, normalized to the maximal activity obtained. The activity assays were performed on two separate days; therefore, data sets for each biological triplicate were scaled according to the activity at 0.35% and then normalized to the maximum activity. Points represent the average of normalizations, and error bars represent standard deviation.



**Figure S2.** The activity of β-gal associated with different aTc plasmid constructs in *H. modesticaldum* with and without aTc inducer (200 µg L<sup>-1</sup>). The constructs include P<sub>gapDH/tetO</sub> controlling lacZ (left), and a second construct including the tetR repressor controlled by P<sub>miniph1</sub> (right). Data represent the average rate of biological and technical triplicate; error bars represent the average standard deviation of each set of biological replicates. The ordinate scale is expanded 100-fold in plot B in order to see the activity of cells harboring the second construct.

## References

1. Olson, D.G.; Maloney, M.; Lanahan, A.A.; Hon, S.; Hauser, L.J.; Lynd, L.R. Identifying Promoters for Gene Expression in *Clostridium Thermocellum*. *Metab. Eng. Commun.* **2015**, *2*, 23–29, doi:10.1016/j.meten.2015.03.002.
2. Nariya, H.; Miyata, S.; Kuwahara, T.; Okabe, A. Development and Characterization of a Xylose-Inducible Gene Expression System for *Clostridium Perfringens*. *J. Microbiol. Methods* **2011**, *77*, 8439–8441, doi:10.1128/AEM.05668-11.
3. Heap, J.T.; Pennington, O.J.; Cartman, S.T.; Minton, N.P. A Modular System for Clostridium Shuttle Plasmids. *J. Microbiol. Methods* **2009**, *78*, 79–85, doi:10.1016/j.mimet.2009.05.004.