

Electrochemical sensing of dopamine using polypyrrole/Molybdenum Oxide modified ITO electrode

Nadiyah Alahmadi ¹ and Waleed Ahmed El-Said ^{1,2 *}

¹ Department of Chemistry, College of Science, University of Jeddah, Jeddah 21959, Saudi Arabia

² Chemistry Department, Faculty of Science, Assiut University, Assiut 71516, Egypt

* Correspondence: awaleedahmed@yahoo.com




Figure S1. Distribution of the particle sizes of (a) MoO₃ NPs/ITO, and (b) PPy/MoO₃ NPs/ITO.