

Research Article

Estimation, evaluation and characterization of carbapenem resistance burden from a tertiary care hospital, Pakistan

Aamir Jamal Gondal¹, Nakhshab Choudhry², Hina Bukhari³, Zainab Rizvi⁴, Shah Jahan⁵ and Nighat Yasmin¹

¹ Department of Biomedical Sciences, King Edward Medical University, Lahore 54000, Pakistan.

² Department of Biochemistry, King Edward Medical University, Lahore 54000, Pakistan.

³ Department of Pathology, King Edward Medical University, Lahore 54000, Pakistan.

⁴ Department of Oral Pathology, de'Montmorency College of Dentistry, Lahore 54000, Pakistan.

⁵ Department of Immunology, University of Health Sciences, Lahore 54600, Pakistan.

Correspondence: Dr. Nighat Yasmin, Associate Professor

Postal Address: Department of Biomedical Sciences, King Edward Medical University, Nelagumbad Chowk, Anarkali Bazar, Lahore 54000, Pakistan.

E.mail: drnyasmin@kemu.edu.pk

Cell: +92-321-8566886

Table S1: Primer sequences for PCR and sequencing analysis

	Primer Sequences	Tm °C / Amplicon (bp)
Antimicrobial Resistance Genes [14, 104, 105]:		
PCR cycling condition: 40 sec at 95 °C, 30 sec at melting temperature, 30 sec at 72 °C		
<i>blaKPC-2</i>	F: gct aca cct agc tcc acc ttc R: aca gtg gtt ggt aat cca tgc	55 / 989
<i>blaNDM-1</i>	F: ggg cag tcg ctt cca acg gt R: gta gtg ctc agt gtc ggc at	53 / 476
<i>blaNDM-1</i>	F: atg gaa ttg ccc aat att atg c R: tca gcg cag ctt gtc ggc cat	50 / 813
<i>blaOXA-48</i>	F: gcg tgg tta agg atg aac ac R: cat caa gtt caa ccc aac cg	52 / 438
<i>blavIM</i>	F: gat ggt gtt tgg tcg cat a R: cga atg cgc agc acc ag	52 / 390
<i>blaIMP</i>	F: gga ata gag tgg ctt aay tct c R: ggt tta aya aaa caa cca cc	52 / 232
<i>blasHV</i>	F: ctt tat cgg ccc tca ctc aa R: agg tgc tca tca tgg gaa ag	55 / 237
<i>blaTEM</i>	F: cgc cgc ata cac tat tct cag aat ga R: acg ctc acc ggc tcc aga ttt at	55 / 445
<i>blaCTX-M</i>	F: atg tgc agy acc agt aar gtk atg gc R: tgg gtr aar tar gts acc aga ayc agc gg	55 / 593
MLST <i>K. pneumoniae</i> [101]: 30 sec at 96 °C, 30 sec at melting temperature, 40 sec at 72 °C		
<i>gapA</i>	F: tga agt atg act cca ctc acg g R: aac gcc ttt cat tgc gcc ttc gga a	60 / 662
<i>infB</i>	F: ctc tct gct gga cta cat tcg R: cgc ttt cag ctc cag aac ttc	52 / 462
<i>mdh</i>	F: ccc aac tgc ctt cag gtt cag R: cct tcc acg tag gcg cat tcc	52 / 756
<i>pgi</i>	F: gag aaa aac ctg cgc gtg ctg ctg R: cgg tta atc agg cgc tta gtg gag c	52 / 566
<i>phoE</i>	F: acc tgg cgc aac acc gat ttc ttc R: ttc agc tgg ttg att ttg taa tcc ac	52 / 602
<i>rpoB</i>	F: ggc gaa atg gcg gaa aac ca R: gag tct tcg aag ttg taa cc	52 / 1075
<i>tonB</i>	F: ctc tat act tcg gta cat cag gtt R: cct gtt tgg cgg cca gca cct ggt	48 / 539
<i>infB2</i>	F: act aag gtt gcc tcc ggc gaa gc	60
<i>pgi2</i>	F: ctg ctg gcg ctg atc ggc at	60
<i>pgi2</i>	R: tta tag cgg tta atc agg cgc t	60
MLST <i>E. coli</i> [102]: PCR cycling condition: 30 sec at 95 °C, 40 sec at 55 °C, 40 sec at 72 °C		
<i>dinB</i>	F: gtt ttc cca gtc acg acg ttg tat gag agg tga gca atg cgt a R: ttg tga gcg gat aac aat ttc cgt agc ccc atc gct tcc ag	
<i>icd2</i>	F: gtt ttc cca gtc acg acg ttg taa ttc gct tcc cgg aac att g R: ttg tga gcg gat aac aat ttc atg atc gcg tca cca aay tc	
<i>pabB</i>	F: gtt ttc cca gtc acg acg ttg taa atc caa tat gac cgg cga g R: ttg tga gcg gat aac aat ttc ggt tcc agt tcg tcc ata at	
<i>polB</i>	F: gtt ttc cca gtc acg acg ttg tag gcg gct atg tga tgg att c R: ttg tga gcg gat aac aat ttc ggt tgg cat cag aaa acg gc	
<i>putP</i>	F: gtt ttc cca gtc acg acg ttg tac tgt tta acc cgt gga ttg c R: ttg tga gcg gat aac aat ttc gca tcg gcc tcg gca aag cg	
<i>trpA</i>	F: gtt ttc cca gtc acg acg ttg tag cta cga atc tct gtt tgc c R: ttg tga gcg gat aac aat ttc gct ttc atc ggt tgt aca aa	
<i>trpB</i>	F: gtt ttc cca gtc acg acg ttg tac act ata tgc tgg gca ccc c	

uidA R: ttg tga gcg gat aac aat ttc cct cgt gct ttc aaa ata tc
F: gtt ttc cca gtc acg acg ttg tac att acg gca aag tgt ggg tca at
R: ttg tga gcg gat aac aat ttc cca tca gca cgt tat cga atc ctt

Sequencing
oF: gtt ttc cca gtc acg acg ttg ta
oR: ttg tga gcg gat aac aat ttc
